



AUTHORIZED DISTRIBUTOR



We've got you covered.

Product Catalog | Acute Care

For those on the front line of patient care.

Our proven, high-quality, single-use infection control solutions are expertly designed for the protection and comfort you deserve, and compliance you can count on.

Lending peace of mind so you can focus on those you care for.



Our Portfolio

Respirators and Masks	3
Gloves	9
Gowns and Apparel	18
Wound Care	28
Sterilization	45
Essential Supplies.....	48



Respirators & Masks

Fit for the task.

Our innovative designs and premium materials ensure that Medicom face masks provide reliable protection, as well as exceptional breathability and comfort, whether NIOSH-certified N95 respirator or ASTM level 1, 2, or 3 face masks.



SafeMask Architect Pro™ Surgical N95 Respirator

Medical grade with optimal respiratory protection.

Filters 95% of particulate aerosols free of oil.

- High-quality materials and duckbill shape offers optimal breathability with no exhalation valve needed.
- Unique single-split comfort bands hold respirator in place without pulling or irritating skin.
- Lightweight construction ensures comfort; up to 50% lighter than other leading respirators. *Based on internal average mass weight assessment vs other respirators on the market.*

Product Code	Description	Package/Case
203214	Small, White	○ 50/box, 10 boxes/case, 500/case
203314	Medium, White	○ 50/box, 10 boxes/case, 500/case
203414	Large, White	○ 50/box, 10 boxes/case, 500/case



AssureMask Balance® Earloop Mask

Product Code	Description	Package/Case
205415	ASTM Level 1, Blue	● 50/box, 10 boxes/case, 500/case
205417	ASTM Level 1, Yellow	● 50/box, 10 boxes/case, 500/case
205416	ASTM Level 1, Pink	● 50/box, 10 boxes/case, 500/case
205414	ASTM Level 1, White	○ 50/box, 10 boxes/case, 500/case
205515	ASTM Level 2, Blue	● 50/box, 10 boxes/case, 500/case
205517	ASTM Level 2, Yellow	● 50/box, 10 boxes/case, 500/case
205516	ASTM Level 2, Pink	● 50/box, 10 boxes/case, 500/case
205514	ASTM Level 2, White	○ 50/box, 10 boxes/case, 500/case
205615	ASTM Level 3, Blue	● 50/box, 10 boxes/case, 500/case
205617	ASTM Level 3, Yellow	● 50/box, 10 boxes/case, 500/case
205616	ASTM Level 3, Pink	● 50/box, 10 boxes/case, 500/case
205614	ASTM Level 3, White	○ 50/box, 10 boxes/case, 500/case



AssureMask Valiant™ Earloop Mask

*Powerful protection with customized comfort.
Made in North America. [Levels 1-2-3]*

- Shingle pleat design prevents fluid pooling and cross-contamination.
- The adjustable, aluminum nose piece forms a strong seal for maximum protection.

Product Code	Description	Package/Case
205715	ASTM Level 1, Blue	50/box, 10 boxes/case, 500/case
205717	ASTM Level 1, Yellow	50/box, 10 boxes/case, 500/case
205716	ASTM Level 1, Pink	50/box, 10 boxes/case, 500/case
205714	ASTM Level 1, White	50/box, 10 boxes/case, 500/case
205815	ASTM Level 2, Blue	50/box, 10 boxes/case, 500/case
205817	ASTM Level 2, Yellow	50/box, 10 boxes/case, 500/case
205816	ASTM Level 2, Pink	50/box, 10 boxes/case, 500/case
205814	ASTM Level 2, White	50/box, 10 boxes/case, 500/case
205915	ASTM Level 3, Blue	50/box, 10 boxes/case, 500/case
205917	ASTM Level 3, Yellow	50/box, 10 boxes/case, 500/case
205916	ASTM Level 3, Pink	50/box, 10 boxes/case, 500/case
205914	ASTM Level 3, White	50/box, 10 boxes/case, 500/case



Ritmed® ProShield with Visor Earloop Anti-fog Masks

Product Code	Description	Package/Case
2125	ASTM Level 3, Teal	25/box, 4 boxes/case, 100/case

Tie-on Anti-fog Masks

2126	ASTM Level 3, Blue	25/box, 4 boxes/case, 100/case
------	--------------------	--------------------------------



AssureMask Precision™ Surgical Tie-On Masks with SureView™ Technology

Fog-free protection with a personalized fit. [Levels 1-2-3]

- A soft, comfortable, and humidity absorbent anti-fog nose strip ensures an unobstructed view of the operative field.
- Single unit dispensing.
- Precision™ Care option includes a soft inner lining made especially for sensitive skin.

Product Code	Description	Package/Case
206815	ASTM Level 1, Blue	50/box, 6 boxes/case, 300/case
206818	ASTM Level 1, Green	50/box, 6 boxes/case, 300/case
206814	ASTM Level 1, White	50/box, 6 boxes/case, 300/case
207015	ASTM Level 2, Blue	50/box, 6 boxes/case, 300/case
207018	ASTM Level 2, Green	50/box, 6 boxes/case, 300/case
207014	ASTM Level 2, White	50/box, 6 boxes/case, 300/case
207115	ASTM Level 3, Blue	50/box, 6 boxes/case, 300/case
207118	ASTM Level 3, Green	50/box, 6 boxes/case, 300/case
207114	ASTM Level 3, White	50/box, 6 boxes/case, 300/case

AssureMask Precision™ Care Surgical Tie-on Mask

206914	ASTM Level 1, White	50/box, 6 boxes/case, 300/case
--------	---------------------	--------------------------------



SafeMask Premier Elite ProShield™ Earloop Mask with Visor

Product Code	Description	Package/Case
204215	ASTM Level 3, Blue	25/box, 4 boxes/case, 100/case
204216	ASTM Level 3, Pink	25/box, 4 boxes/case, 100/case



SafeMask TailorMade™ Earloop Mask

Product Code	Description	Package/Case
2070	ASTM Level 1, Blue	50/box, 10 boxes/case, 500/case
2071	ASTM Level 1, Pink	50/box, 10 boxes/case, 500/case
2072	ASTM Level 3, Blue	50/box, 10 boxes/case, 500/case



SafeMask FreeFlow® with Air Space™ Technology

Omega pleat design creates best-in-class breathability and comfort. [Levels 1-2-3]

- Innovative contoured pleats create a unique pocket that allows air to flow freely.
- Integrated strip ensures clear vision preventing eyewear from fogging.

Product Code	Description	Package/Case
200311	ASTM Level 1, Lavender	50/box, 10 boxes/case, 500/case
200315	ASTM Level 1, Blue	50/box, 10 boxes/case, 500/case
200316	ASTM Level 1, Pink	50/box, 10 boxes/case, 500/case
200415	ASTM Level 2, Blue	50/box, 10 boxes/case, 500/case
200416	ASTM Level 2, Pink	50/box, 10 boxes/case, 500/case
200511	ASTM Level 3, Lavender	50/box, 10 boxes/case, 500/case
200514	ASTM Level 3, White	50/box, 10 boxes/case, 500/case
200515	ASTM Level 3, Blue	50/box, 10 boxes/case, 500/case
200516	ASTM Level 3, Pink	50/box, 10 boxes/case, 500/case



Midnight™ Earloop Mask

Product Code	Description	Package/Case
213792	ASTM Level 1, Black	50/box, 10 boxes/case, 500/case
213892	ASTM Level 2, Black	50/box, 10 boxes/case, 500/case
213992	ASTM Level 3, Black	50/box, 10 boxes/case, 500/case



SafeMask® Junior Earloop Pediatric Mask

Product Code	Description	Package/Case
2035	ASTM Level 1, Blue	50/box, 10 boxes/case, 500/case



SafeMask SofSkin™ Earloop Mask

*Reliable protection made especially for sensitive skin.
[Levels 1 and 3]*

- Lightweight, breathable inner layer softens with wear.
- Unique wicking technology keeps face dry and comfortable.
- Ideal for those with skin sensitivity.

Product Code	Description	Package/Case
2080	ASTM Level 1, White	○ 50/box, 10 boxes/case, 500/case
2081	ASTM Level 1, Blue	● 50/box, 10 boxes/case, 500/case
2082	ASTM Level 1, Pink	● 50/box, 10 boxes/case, 500/case
2083	ASTM Level 1, Lavender	● 50/box, 10 boxes/case, 500/case
2086	ASTM Level 3, White	○ 50/box, 10 boxes/case, 500/case
2087	ASTM Level 3, Blue	● 50/box, 10 boxes/case, 500/case
2088	ASTM Level 3, Pink	● 50/box, 10 boxes/case, 500/case
2089	ASTM Level 3, Lavender	● 50/box, 10 boxes/case, 500/case



SafeMask® Cone Mask

Product Code	Description	Package/Case
2020	Adjustable Aluminium Nose Piece, Blue	● 50/box, 12 boxes/case, 600/case



Mask Box Holder

Product Code	Description	Package/Case
MAS001	Box holder / dispenser masks	2/pack



Gloves

A perfect pair.

Our complete range of premium quality medical examination and surgical gloves are made using select materials that provide superior protection, durability, and grip without compromise to flexibility, tactility, or comfort.



AssureTouch Agile™

Nitrile glove for mobility and protection.

- Superior puncture, tear, and chemical resistance without compromising dexterity or tactility.
- Textured fingertips ensure a secure grip, even when wet.
- Chemo tested.



Nitrile Glove
Chemo Tested, Powder-free, 4.5 mil

Product Code	Description	Package/Case
9998-A	XS, Blue	● 100/box, 10 boxes/case, 1000/case
9998-B	S, Blue	● 100/box, 10 boxes/case, 1000/case
9998-C	M, Blue	● 100/box, 10 boxes/case, 1000/case
9998-D	L, Blue	● 100/box, 10 boxes/case, 1000/case
9998-E	XL, Blue	● 100/box, 10 boxes/case, 1000/case



AssureTouch™ Pure Nitrile Glove
DermaSense Technology, Chemo Tested, Powder-free, 2.76 mil

Product Code	Description	Package/Case
9970-A	XS, Indigo	● 200/box, 10 boxes/case, 2000/case
9970-B	S, Indigo	● 200/box, 10 boxes/case, 2000/case
9970-C	M, Indigo	● 200/box, 10 boxes/case, 2000/case
9970-D	L, Indigo	● 200/box, 10 boxes/case, 2000/case
9970-E	XL, Indigo	● 200/box, 10 boxes/case, 2000/case



AssureTouch Ultra Tough™ Nitrile Glove
Fully Textured, Powder-free, 5 mil

Product Code	Description	Package/Case
9950-B	S, Black	● 100/box, 10 boxes/case, 1000/case
9950-C	M, Black	● 100/box, 10 boxes/case, 1000/case
9950-D	L, Black	● 100/box, 10 boxes/case, 1000/case
9950-E	XL, Black	● 100/box, 10 boxes/case, 1000/case
9950-F	XXL, Black	● 100/box, 10 boxes/case, 1000/case



AssureTouch™ Protect

Nitrile glove for reliable protection and enhanced sensitivity.

- Thin nitrile glove ideal for procedures that require precision and tactile sensitivity.
- Chemo tested.



Nitrile Glove
Chemo Tested, Powder-free, 3 mil

Product Code	Description	Package/Case
9990-A	XS, Blue	● 100/box, 10 boxes/case, 1000/case
9990-B	S, Blue	● 100/box, 10 boxes/case, 1000/case
9990-C	M, Blue	● 100/box, 10 boxes/case, 1000/case
9990-D	L, Blue	● 100/box, 10 boxes/case, 1000/case
9990-E	XL, Blue	● 100/box, 10 boxes/case, 1000/case



AssureTouch Power™ Nitrile Glove
Extended Cuff, Fentanyl Safe, Chemo Tested, Powder-free, 4 mil

Product Code	Description	Package/Case
9999-A	XS, Blue	● 100/box, 10 boxes/case, 1000/case
9999-B	S, Blue	● 100/box, 10 boxes/case, 1000/case
9999-C	M, Blue	● 100/box, 10 boxes/case, 1000/case
9999-D	L, Blue	● 100/box, 10 boxes/case, 1000/case
9999-E	XL, Blue	● 100/box, 10 boxes/case, 1000/case



AssureTouch™ Basics Nitrile Glove
Chemo Tested, Powder-free, 2.7 mil

Product Code	Description	Package/Case
9985-A	XS, Blue	● 100/box, 10 boxes/case, 1000/case
9985-B	S, Blue	● 100/box, 10 boxes/case, 1000/case
9985-C	M, Blue	● 100/box, 10 boxes/case, 1000/case
9985-D	L, Blue	● 100/box, 10 boxes/case, 1000/case
9985-E	XL, Blue	● 100/box, 10 boxes/case, 1000/case



AssureTouch Respond™ Latex Glove

Fully Textured, Powder-free, 4.7 mil

Product Code	Description	Package/Case
9992-A	XS, Natural	● 100/box, 10 boxes/case, 1000/case
9992-B	S, Natural	● 100/box, 10 boxes/case, 1000/case
9992-C	M, Natural	● 100/box, 10 boxes/case, 1000/case
9992-D	L, Natural	● 100/box, 10 boxes/case, 1000/case
9992-E	XL, Natural	● 100/box, 10 boxes/case, 1000/case



AssureTouch Shield™ Vinyl Glove

3.2 mil

Product Code	Description	Package/Case
9994-A	S, Clear	○ 100/box, 10 boxes/case, 1000/case
9994-B	M, Clear	○ 100/box, 10 boxes/case, 1000/case
9994-C	L, Clear	○ 100/box, 10 boxes/case, 1000/case
9994-D	XL, Clear	○ 100/box, 10 boxes/case, 1000/case



Ritmed® Vital Touch™ Latex Procedure Glove

Sterile, Powder-free, 3.9 mil

Product Code	Description	Package/Case
1143-6	6, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-6.5	6.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-7	7, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-7.5	7.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-8	8, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-8.5	8.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1143-9	9, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case



Ritmed® Vital Touch™

Sterile premium latex surgical glove.

- Optimal dexterity and tactility while providing reliable barrier protection.
- Extended cuffs allow quick and easy donning.
- Each pair is individually sealed to maintain sterility.



Latex Surgical Glove Sterile, Powder-free, 4 mil

Product Code	Description	Package/Case
1141-5.5	5.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-6	6, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-6.5	6.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-7	7, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-7.5	7.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-8	8, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case
1141-8.5	8.5, White	○ 1 pair/pouch, 50 pouches/box, 8 boxes/case, 400 pairs/case



SafeTouch MicroDefense™ Antimicrobial Nitrile Glove Chemo Tested, Powder-free, 2.5 mil

Product Code	Description	Package/Case
1230-A	XS, Violet Blue	● 200/box, 10 boxes/case, 2000/case
1230-B	S, Violet Blue	● 200/box, 10 boxes/case, 2000/case
1230-C	M, Violet Blue	● 200/box, 10 boxes/case, 2000/case
1230-D	L, Violet Blue	● 200/box, 10 boxes/case, 2000/case
1230-E	XL, Violet Blue	● 180/box, 10 boxes/case, 1800/case



SafeTouch Rejuvenate™ Latex Glove
Lanolin/Vitamin E, Powder-free, 5.1 mil

Product Code	Description	Package/Case
1163-A	XS, Natural	100/box, 10 boxes/case, 1000/case
1163-B	S, Natural	100/box, 10 boxes/case, 1000/case
1163-C	M, Natural	100/box, 10 boxes/case, 1000/case
1163-D	L, Natural	100/box, 10 boxes/case, 1000/case
1163-E	XL, Natural	100/box, 10 boxes/case, 1000/case



SafeBasics True Fit™ Thin Nitrile Glove
Fentanyl Safe, Chemo Tested, Powder-free, 3 mil

Product Code	Description	Package/Case
1185-A	XS, Blue	300/box, 8 boxes/case, 2400/case
1185-B	S, Blue	300/box, 8 boxes/case, 2400/case
1185-C	M, Blue	300/box, 8 boxes/case, 2400/case
1185-D	L, Blue	300/box, 8 boxes/case, 2400/case
1185-E	XL, Blue	300/box, 8 boxes/case, 2400/case

Chemo Tested, Powder-free, 2.5 mil

1186-A	XS, Pink	300/box, 8 boxes/case, 2400/case
1186-B	S, Pink	300/box, 8 boxes/case, 2400/case
1186-C	M, Pink	300/box, 8 boxes/case, 2400/case
1186-D	L, Pink	300/box, 8 boxes/case, 2400/case



Ocean Pacific® Intuition® Nitrile Glove
Chemo Tested, Powder-free, 2.4 mil

Product Code	Description	Package/Case
OPXS-IN	XS, Lilac	200/box, 10 boxes/case, 2000/case
OPS-IN	S, Lilac	200/box, 10 boxes/case, 2000/case
OPM-IN	M, Lilac	200/box, 10 boxes/case, 2000/case
OPL-IN	L, Lilac	200/box, 10 boxes/case, 2000/case
OPXL-IN	XL, Lilac	200/box, 10 boxes/case, 2000/case



Ocean Pacific® Lab Master®

DermaSense Technology helps reduce the risk of skin reactions.

- Protection and strength of nitrile without chemical accelerators.
- Optimal tactile sensitivity.
- Chemo tested.



Nitrile Glove
Chemo Tested, Powder-free, 2.76 mil

Product Code	Description	Package/Case
OPXS-LM2	XS, Indigo	● 200/box, 10 boxes/case, 2000/case
OPS-LM2	S, Indigo	● 200/box, 10 boxes/case, 2000/case
OPM-LM2	M, Indigo	● 200/box, 10 boxes/case, 2000/case
OPL-LM2	L, Indigo	● 200/box, 10 boxes/case, 2000/case
OPXL-LM2	XL, Indigo	● 200/box, 10 boxes/case, 2000/case



Ocean Pacific® Blue Lite® Nitrile Glove
Chemo Tested, Powder-free, 2.8 mil

Product Code	Description	Package/Case
OPXS-BL2	XS, Blue	● 200/box, 10 boxes/case, 2000/case
OPS-BL2	S, Blue	● 200/box, 10 boxes/case, 2000/case
OPM-BL2	M, Blue	● 200/box, 10 boxes/case, 2000/case
OPL-BL2	L, Blue	● 200/box, 10 boxes/case, 2000/case
OPXL-BL2	XL, Blue	● 200/box, 10 boxes/case, 2000/case



Ocean Pacific® NeoNatural™ Polychloroprene Glove
Chemical Resistant, Chemo Tested, 4 mil

Product Code	Description	Package/Case
1228-A	XS, Emerald	● 100/box, 10 boxes/case, 1000/case
1228-B	S, Emerald	● 100/box, 10 boxes/case, 1000/case
1228-C	M, Emerald	● 100/box, 10 boxes/case, 1000/case
1228-D	L, Emerald	● 100/box, 10 boxes/case, 1000/case
1228-E	XL, Emerald	● 100/box, 10 boxes/case, 1000/case



Ocean Pacific® Eclipse® Nitrile Glove
Powder-free, 1.96 mil

Product Code	Description	Package/Case
OPXS-EC	XS, Blue	● 200/box, 10 boxes/case, 2000/case
OPS-EC	S, Blue	● 200/box, 10 boxes/case, 2000/case
OPM-EC	M, Blue	● 200/box, 10 boxes/case, 2000/case
OPL-EC	L, Blue	● 200/box, 10 boxes/case, 2000/case
OPXL-EC	XL, Blue	● 200/box, 10 boxes/case, 2000/case



Ocean Pacific® Ninja® Nitrile Glove
Fentanyl Safe, Chemo Tested, Fully Textured, Powder-free, 4.8 mil

Product Code	Description	Package/Case
OPS-NI	S, Black	● 200/box, 10 boxes/case, 2000/case
OPM-NI	M, Black	● 200/box, 10 boxes/case, 2000/case
OPL-NI	L, Black	● 200/box, 10 boxes/case, 2000/case
OPXL-NI	XL, Black	● 200/box, 10 boxes/case, 2000/case



General Purpose Glove Options
Octo™ Copolymer, Textured, Powder-free, 4.3 mil

Product Code	Description	Package/Case
9925-B	S, Blue	● 100/box, 10 boxes/case, 1000/case
9925-C	M, Blue	● 100/box, 10 boxes/case, 1000/case
9925-D	L, Blue	● 100/box, 10 boxes/case, 1000/case
9925-E	XL, Blue	● 100/box, 10 boxes/case, 1000/case
9925-F	XXL, Blue	● 100/box, 10 boxes/case, 1000/case
9935-B	S, Clear	○ 100/box, 10 boxes/case, 1000/case
9935-C	M, Clear	○ 100/box, 10 boxes/case, 1000/case
9935-D	L, Clear	○ 100/box, 10 boxes/case, 1000/case
9935-E	XL, Clear	○ 100/box, 10 boxes/case, 1000/case
9935-F	XXL, Clear	○ 100/box, 10 boxes/case, 1000/case

Polyethylene, Embossed, Powder-free, 0.3 mil

9945-B	S, Clear	○ 500/box, 20 boxes/case, 10000/case
9945-C	M, Clear	○ 500/box, 20 boxes/case, 10000/case
9945-D	L, Clear	○ 500/box, 20 boxes/case, 10000/case
9945-E	XL, Clear	○ 500/box, 20 boxes/case, 10000/case



Long Sleeved Insemination Glove 1.18 mil

Product Code	Description	Package/Case
9965	Universal Size, Red	● 150/box, 10 boxes/case, 1500/case
9966	Universal Size, Yellow	● 100/box, 10 boxes/case, 1000/case



Glove Box Holders

Product Code	Description	Package/Case
GLV001	Single Glove Box Holder	2/case
GLV003	Triple Glove Box Holder	1/case



Gowns & Apparel

We've got your back.

Our complete line of medical gowns and protective apparel are made of premium materials and designed to provide reliable barrier protection with the added benefit of comfort and style.



AssureWear™ Sterile Surgical Gown

Sterile ultrasonically sealed gown for reliable protection.

- Latex-free fabric and adhesive neck closure minimize the risk of skin reactions.
- Removable by pulling forward at shoulder level.
- Two 4-ply absorbent towels included.



AAMI Level 3

Multi-layer SMS, Velcro Neck, Knit Cuff, Two Ties Back Closure

Product Code	Description	Package/Case
S67000	M, Blue	1/bag, 50 bags/case
S67001	L, Blue	1/bag, 50 bags/case
S67002	XL, Blue	1/bag, 50 bags/case



AAMI Level 3 - AssureWear™ Chemo-Tested VersaGown® with Flexneck™ Technology

Non-sterile, Poly Coated, Thumb Loop, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69952	XL, Yellow	10/bag, 10 bags/case, 100/case
A69956	XL, Blue	10/bag, 10 bags/case, 100/case



AAMI Level 2 - AssureWear™ Chemo-Tested VersaGown® with Flexneck Plus™ Technology

Non-sterile, Poly Coated, Knit Cuffs, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69931	L, Blue	10/bag, 10 bags/case, 100/case
A69932	XL, Blue	10/bag, 10 bags/case, 100/case



AssureWear™ Chemo-Tested VersaGown® with Flexneck™ Technology

Reduces cross-contamination and increases compliance.

- Double elastic neck closure allows quick and safe donning and doffing.
- Thumb loop cuffs prevent sleeves from riding up.
- Clinically proven to reduce the risk of cross-contamination versus standard gowns. *Data on file.*



AAMI Level 2 Isolation Gown Non-sterile, Poly Coated, Thumb Loop, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69949	S, Blue	10/bag, 10 bags/case, 100/case
A69950	M, Blue	10/bag, 10 bags/case, 100/case
A69954	L, Blue	10/bag, 10 bags/case, 100/case
A69955	XL, Blue	10/bag, 10 bags/case, 100/case
A69951	L, Yellow	10/bag, 10 bags/case, 100/case



AAMI Level 2 - AssureWear™ Chemo-Tested Non-sterile, Poly Coated, Neck Ties, Thumb Loop, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69954-T	L, Blue	10/bag, 10 bags/case, 100/case



AAMI Level 2 - AssureWear™ VersaGown® with Flexneck™ Technology Non-sterile, Multi-layer SMS, Thumb Loop, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69921	L, Yellow	10/bag, 10 bags/case, 100/case
A69922	XL, Yellow	10/bag, 10 bags/case, 100/case
A69961	L, Yellow	10/bag, 10 bags/case, 100/case
A69962	XL, Yellow	10/bag, 10 bags/case, 100/case
A69963	2XL, Yellow	10/bag, 10 bags/case, 100/case



Ritmed® VersaGown Agile™ with Flexneck Plus™ and TieFree™ Technologies

Promotes instant infection control compliance.

- Double elastic neck closure, overlapping back and tie-free design saves time when donning and doffing.
- Thumb loop cuffs prevent sleeves from riding up.
- Clinically proven to reduce the risk of cross-contamination versus standard gowns. *Data on file.*



AAMI Level 2 Isolation Gown Non-sterile, Multi-layer SMS, Thumb Loop, Full Back

Product Code	Description	Package/Case
A69933	L, Yellow	10/bag, 10 bags/case, 100/case
A69934	XL, Yellow	10/bag, 10 bags/case, 100/case



AAMI Level 2 - AssureWear™ VersaGown® with Flexneck™ Technology Non-sterile, Multi-layer SMS, Elastic Cuffs, Full Back, Belt on the Front

Product Code	Description	Package/Case
A69980	S, White	10/bag, 10 bags/case, 100/case
A69981	M, White	10/bag, 10 bags/case, 100/case
A69982	XL, White	10/bag, 10 bags/case, 100/case
A69983	3XL, White	10/bag, 10 bags/case, 100/case
A69984	7XL, White	10/bag, 10 bags/case, 100/case

Multi-Layer SMS, Thumb Loop, Full Back, Belt on the Front

A69979-OTH	L, Yellow	10/bag, 10 bags/case, 100/case
------------	-----------	--------------------------------



AAMI Level 2 - AssureWear™ Isolation Gown Non-sterile, Multi-layer SMS, Neck Ties, Elastic Cuffs, Full Back, Belt on the Front

Product Code	Description	Package/Case
A69972	L, Yellow	10/bag, 10 bags/case, 100/case
A69973	XL, Yellow	10/bag, 10 bags/case, 100/case



AssureWear™ VersaGown® with Flexneck™ Technology

Extended thumb loop cuffs further reduce the risk of cross-contamination.

- Double elastic neck closure allows quick and safe donning and doffing.
- Extended thumb loop cuffs increase coverage of the palms and eliminate gaping at the wrists.
- Clinically proven to reduce the risk of cross-contamination versus standard gowns. *Data on file.*



AAMI Level 2 - Isolation Gown Non-sterile, Multi-layer SMS, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69927	L, Yellow	10/bag, 10 bags/case, 100/case
A69928	XL, Yellow	10/bag, 10 bags/case, 100/case



AAMI Level 2 - AssureWear™ Isolation Gown Non-sterile, Multi-layer SMS, Over-The-Head, Thumb Loop, Full Back, Belt on the Left

Product Code	Description	Package/Case
A69977	L, Yellow	10/bag, 10 bags/case, 100/case
A69978	XL, Yellow	10/bag, 10 bags/case, 100/case



AAMI Level 1 - AssureWear™ VersaGown® with Flexneck™ Technology Non-sterile, Multi-layer SMS, Thumb Loop, Full Back, Belt on the Front

Product Code	Description	Package/Case
A69911	L, Yellow	10/bag, 10 bags/case, 100/case
A69912	XL, Yellow	10/bag, 10 bags/case, 100/case
A69915	L, Blue	10/bag, 10 bags/case, 100/case
A69916	XL, Blue	10/bag, 10 bags/case, 100/case



AAMI Level 1 - Isolation Gown

Non-sterile, Multi-layer SMS, Neck Ties, Thumb Loop, Full Back, Belt on the Front

Product Code	Description	Package/Case
A69911-T	L, Yellow	10/bag, 10 bags/case, 100/case
A69912-T	XL, Yellow	10/bag, 10 bags/case, 100/case




Non-AAMI Rated - Ritmed® Isolation Gown

Non-sterile, Spunbond Polypropylene, Neck Ties, Elastic Cuffs, Full Back, Belt on the Right

Product Code	Description	Package/Case
D-8009	Reg, Yellow	10/bag, 5 bags/case, 50/case
D-8010	XL, Yellow	10/bag, 5 bags/case, 50/case

Polyethylene, Over-the-Head, Thumb Loop, Open Back


8025	Reg, Blue	1/bag, 15 bags/box, 5 boxes/case, 75/case
8026	XL, Blue	1/bag, 15 bags/box, 5 boxes/case, 75/case



Non-AAMI Rated - SafeWear™ Isolation Gown

Non-sterile, Spunbond Polypropylene, Neck Ties, Knit Cuffs, Full Back, Belt on the Right

Product Code	Description	Package/Case
8113	Reg, Sunny Yellow	12/bag, 5 bags/case, 60/case
8113-D	XL, Sunny Yellow	12/bag, 5 bags/case, 60/case
8114	Reg, Bright Blue	12/bag, 5 bags/case, 60/case
8114-D	XL, Bright Blue	12/bag, 5 bags/case, 60/case



Non-AAMI Rated - SafeBasics™ Isolation Gown

Non-sterile, Poly Coated, Neck Ties, Knit Cuffs, Full Back, Belt on the Right

Product Code	Description	Package/Case
8022	Reg, White	10/bag, 5 bags/case, 50/case
8024	XL, White	10/bag, 5 bags/case, 50/case

Poly Coated, Neck Ties, Elastic Cuffs, Full Back, Belt on the Right

8023	Reg, Blue	10/bag, 5 bags/case, 50/case
------	-----------	------------------------------



Non-AAMI Rated - AssureWear™ Chemo-Tested VersaGown® with Flexneck™ Technology
 Non-sterile, Poly Coated, Thumb Loop, Open Back, Belt on the Left

Product Code	Description	Package/Case
A69959	XL, Blue	10/bag, 10 bags/case, 100/case



SafeBasics™ Lab Gowns
 SMS Material, Knit Collar and Cuffs, Front Snap Button Closure, Two Front Pockets

Product Code	Description	Package/Case
8016	S, Blue	10/bag, 5 bags/case, 50/case
8017	M, Blue	10/bag, 5 bags/case, 50/case
8018	L, Blue	10/bag, 5 bags/case, 50/case
8019	XL, Blue	10/bag, 5 bags/case, 50/case



Ritmed® Lab Coats
 SMS Material, Knit Collar and Cuffs, Front Snap Button Closure

Product Code	Description	Package/Case
A8031	S, Blue	10/bag, 5 bags/case, 50/case
A8032	M, Blue	10/bag, 5 bags/case, 50/case
A8033	L, Blue	10/bag, 5 bags/case, 50/case
A8034	XL, Blue	10/bag, 5 bags/case, 50/case
A8035	2XL, Blue	10/bag, 5 bags/case, 50/case
A8036	3XL, Blue	10/bag, 5 bags/case, 50/case
A8038	4XL, Blue	10/bag, 5 bags/case, 50/case
A8039	7XL, Blue	10/bag, 5 bags/case, 50/case
A8040	9XL, Blue	10/bag, 5 bags/case, 50/case
A8041	S, White	10/bag, 5 bags/case, 50/case
A8042	M, White	10/bag, 5 bags/case, 50/case
A8043	L, White	10/bag, 5 bags/case, 50/case
A8044	XL, White	10/bag, 5 bags/case, 50/case
A8045	2XL, White	10/bag, 5 bags/case, 50/case
A8046	3XL, White	10/bag, 5 bags/case, 50/case
A8048	4XL, White	10/bag, 5 bags/case, 50/case
A8049	7XL, White	10/bag, 5 bags/case, 50/case
A8050	9XL, White	10/bag, 5 bags/case, 50/case
A8051	12XL, White	10/bag, 5 bags/case, 50/case



SafeWear™ High Performance Lab Coat
SMS Material, Knit Collar and Cuffs, Front Snap Button Closure, Three Front Pockets

Product Code	Description	Package/Case
8107-A	S, Deep Blue	12/bag, 5 bags/case, 60/case
8107-B	M, Deep Blue	12/bag, 5 bags/case, 60/case
8107-C	L, Deep Blue	12/bag, 5 bags/case, 60/case
8107-D	XL, Deep Blue	12/bag, 5 bags/case, 60/case
8107-E	2XL, Deep Blue	12/bag, 5 bags/case, 60/case
8107-F	3XL, Deep Blue	12/bag, 5 bags/case, 60/case
8108-A	S, Pretty Pink	12/bag, 5 bags/case, 60/case
8108-B	M, Pretty Pink	12/bag, 5 bags/case, 60/case
8108-C	L, Pretty Pink	12/bag, 5 bags/case, 60/case
8108-D	XL, Pretty Pink	12/bag, 5 bags/case, 60/case
8109-A	S, Plum Purple	12/bag, 5 bags/case, 60/case
8109-B	M, Plum Purple	12/bag, 5 bags/case, 60/case
8109-C	L, Plum Purple	12/bag, 5 bags/case, 60/case
8109-D	XL, Plum Purple	12/bag, 5 bags/case, 60/case
8110-A	S, White Frost	12/bag, 5 bags/case, 60/case
8110-B	M, White Frost	12/bag, 5 bags/case, 60/case
8110-C	L, White Frost	12/bag, 5 bags/case, 60/case
8110-D	XL, White Frost	12/bag, 5 bags/case, 60/case
8112-A	S, Soft Blue	12/bag, 5 bags/case, 60/case
8112-B	M, Soft Blue	12/bag, 5 bags/case, 60/case
8112-C	L, Soft Blue	12/bag, 5 bags/case, 60/case
8112-D	XL, Soft Blue	12/bag, 5 bags/case, 60/case



SafeWear™ Hipster Jacket™
SMS Material, Knit Collar and Cuffs, Front Snap Button Closure, Three Front Pockets

Product Code	Description	Package/Case
8101-A	S, Deep Blue	12/bag, 5 bags/case, 60/case
8101-B	M, Deep Blue	12/bag, 5 bags/case, 60/case
8101-C	L, Deep Blue	12/bag, 5 bags/case, 60/case
8101-D	XL, Deep Blue	12/bag, 5 bags/case, 60/case
8101-E	2XL, Deep Blue	12/bag, 5 bags/case, 60/case
8101-F	3XL, Deep Blue	12/bag, 5 bags/case, 60/case

Continues on Page 26 >



SafeWear™ Hipster Jacket™

SMS Material, Knit Collar and Cuffs, Front Snap Button Closure, Three Front Pockets

Product Code	Description	Package/Case
8102-A	S, Pretty Pink	12/bag, 5 bags/case, 60/case
8102-B	M, Pretty Pink	12/bag, 5 bags/case, 60/case
8102-C	L, Pretty Pink	12/bag, 5 bags/case, 60/case
8102-D	XL, Pretty Pink	12/bag, 5 bags/case, 60/case
8103-A	S, Plum Purple	12/bag, 5 bags/case, 60/case
8103-B	M, Plum Purple	12/bag, 5 bags/case, 60/case
8103-C	L, Plum Purple	12/bag, 5 bags/case, 60/case
8103-D	XL, Plum Purple	12/bag, 5 bags/case, 60/case
8104-A	S, Soft Blue	12/bag, 5 bags/case, 60/case
8104-B	M, Soft Blue	12/bag, 5 bags/case, 60/case
8104-C	L, Soft Blue	12/bag, 5 bags/case, 60/case
8104-D	XL, Soft Blue	12/bag, 5 bags/case, 60/case
8105-A	S, White Frost	12/bag, 5 bags/case, 60/case
8105-B	M, White Frost	12/bag, 5 bags/case, 60/case
8105-C	L, White Frost	12/bag, 5 bags/case, 60/case
8105-D	XL, White Frost	12/bag, 5 bags/case, 60/case



SafeBasics™ Bouffant Caps

Spunbonded Polypropylene

Product Code	Description	Package/Case
8000	21", White	100/box, 5 boxes/case, 500/case
8001	21", Blue	100/box, 5 boxes/case, 500/case
8003	24", White	100/box, 5 boxes/case, 500/case
8004	24", Blue	100/box, 5 boxes/case, 500/case



Ritmed® Bouffant Caps

Spunbonded Polypropylene

Product Code	Description	Package/Case
80027-18	27", Blue	100/box, 5 boxes/case, 500/case
8004-18	24", Blue (in Bag)	100/box, 5 boxes/case, 500/case
A8005	24", Green	100/box, 5 boxes/case, 500/case



SafeBasics™ Shoe Covers
Spunbonded Polypropylene, Non-conductive, Non-skid

Product Code	Description	Package/Case
8006	XL, Blue	100/box, 3 boxes/case, 300/case
8007	Reg, Blue	100/box, 3 boxes/case, 300/case
8006XL	XL, Blue	100/box, 3 boxes/case, 300/case

Spunbonded Polypropylene, Non-conductive

8008	Reg, Blue	100/box, 3 boxes/case, 300/case
8008-XL	XL Blue	100/box, 3 boxes/case, 300/case



Ritmed® Shoe Covers
Plastic, Non-skid

Product Code	Description	Package/Case
A8008-PLAS	Reg, Blue	100/box, 5 boxes/case, 500/case



Ritmed® High Top Boot Covers
Poly Coated

Product Code	Description	Package/Case
A8007-HT-XL	XL, Blue	50/box, 3 boxes/case, 150/case



SafeBasics™ Beard Covers
Spunbonded Polypropylene

Product Code	Description	Package/Case
8035	18", White	100/bag, 10 bags/case, 1000/case



Wound Care

Promoting healing.

Our complete portfolio of high-quality wound care products, including gauze, bandages, sponges, and pads, help keep patients safe and comfortable while they heal.



Ritmed® Vital Gauze Premium Non-Woven Sponges

Blended gauze providing excellent absorbency.

- Protects wounds and helps prevent further injury.
- Fits all parts of the body.



Non-sterile

Product Code	Description	Package/Case
A7632	4-ply, 2" x 2"	200/sleeve, 40 sleeves/case, 8000/case
A7633	4-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
A7634	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case

Sterile

A2338	4-ply, 2" x 2"	2/pouch, 25 pouches/tray, 60 trays/case, 3000/case
A2337	4-ply, 4" x 4"	2/pouch, 25 pouches/tray, 24 trays/case, 1200/case



SafeBasics™ Non-Woven Sponges Non-sterile

Product Code	Description	Package/Case
2101-CH	4-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
2103-CH	4-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
2100-CH	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case



SafeGauze® Premium Non-Woven Sponges with Diamond Aperture™ Non-sterile, Made in USA

Product Code	Description	Package/Case
4522	4-ply, 2" x 2"	200/sleeve, 20 sleeves/case, 4000/case
4544	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case


Ritmed® Non-Woven Sponges
 Non-sterile

Product Code	Description	Package/Case
A2101-CH	4-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
A2103-CH	4-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
A2100-CH	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case
A2106-CH	4-ply, 8" x 4"	200/sleeve, 10 sleeves/case, 2000/case

Sterile

A2240	4-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
A3340	4-ply, 3" x 3"	2/pouch, 50 pouches/tray, 24 trays/case, 2400/case
A4440	4-ply, 4" x 4"	2/pouch, 50 pouches/tray, 12 trays/case, 1200/case
A4440-10	4-ply, 4" x 4"	10/pouch, 100 pouches/case, 1000/case
A8840	4-ply, 8" x 4"	2/pouch, 25 pouches/tray, 20 trays/case, 1000/case
A8840-10	4-ply, 8" x 4"	10/pouch, 100 pouches/case, 1000/case
A9220	6-ply, 2" x 2" (I.V.)	2/pouch, 35 pouches/tray, 20 trays/case, 1400/case
A9220-NP	6-ply, 2" x 2" (I.V. Non Slit)	2/pouch, 35 pouches/tray, 20 trays/case, 1400/case
A9210	6-ply, 4" x 4" (Drain)	2/pouch, 25 pouches/tray, 12 trays/case, 600/case
A9210-NP	6-ply, 4" x 4" (Drain Non Slit)	2/pouch, 25 pouches/tray, 12 trays/case, 600/case


General Purpose Non-Woven Sponges
 Non-sterile

Product Code	Description	Package/Case
2102	4-ply, 2" x 2"	200/sleeve, 20 sleeves/case, 4000/case
2104	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case



Ritmed® Premium Cotton Gauze Sponges

Sterile 100% cotton gauze conveniently packaged.

- Soft and highly absorbent sponges ensure patient comfort.
- Convenient individual peel-open packaging ensures safety.



Non-sterile

Product Code	Description	Package/Case
A3001	8-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
A3001-10	8-ply, 2" x 2" (10s)	200/sleeve, 25 sleeves/case, 5000/case
A3003	8-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
A3005	8-ply, 4" x 4"	200/sleeve, 20 sleeves/case, 4000/case
A3002	12-ply, 2" x 2"	200/sleeves, 40 sleeves/case, 8000/case
A3002-10	12-ply, 2" x 2" (10s)	200/sleeve, 40 sleeves/case, 8000/case
A3004	12-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
A3004-10	12-ply, 3" x 3" (10s)	200/sleeve, 20 sleeves/case, 4000/case
A3006	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case
A3010	12-ply, 4" x 4" (Traech)	200/sleeve, 10 sleeves/case, 2000/case
A3009	12-ply, 8" x 4"	200/sleeve, 10 sleeves/case, 2000/case
A3009-10	12-ply, 8" x 4" (10s)	200/sleeve, 10 sleeves/case, 2000/case

Sterile

A2800	8-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
A4800	8-ply, 4" x 4"	2/pouch, 50 pouches/tray, 12 trays/case, 1200/case
A2120	12-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
A3120	12-ply, 3" x 3"	2/pouch, 40 pouches/tray, 30 trays/case, 2400/case
A4412	12-ply, 4" x 4"	1/pouch, 100 pouches/tray, 12 trays/box, 1200/case
A4120	12-ply, 4" x 4"	2/pouch, 25 pouches/tray, 24 trays/case, 1200/case
A4121	12-ply, 4" x 4"	10/pouch, 10 pouches/tray, 12 trays/box, 1200/case
A8120	12-ply, 8" x 4"	2/pouch, 25 pouches/tray, 20 trays/case, 1000/case
A8121	12-ply, 8" x 4"	10/pouch, 10 pouches/tray, 12 trays/box, 1200/case
A8125	12-ply, 8" x 4"	5/pouch, 10 pouches/tray, 20 trays/case, 1000/case
A4161	16-ply, 4" x 4"	10/pouch, 10 pouches/tray, 12 trays/case, 1200/case


Ritmed® All Cotton Gauze Sponges
 Non-sterile

Product Code	Description	Package/Case
310030	12-ply, 4" x 4"	10 sponges/group, 20 groups/bag, 200 sponges/bag, 10 bags/case, 2000 sponges/case
310040	16-ply, 4" x 4"	10 sponges/group, 20 groups/bag, 200 sponges/bag, 10 bags/case, 2000 sponges/case


Ritmed® Cotton Gauze Sponges
 Non-sterile

Product Code	Description	Package/Case
B3001	8-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
B3001-2	8-ply, 2" x 2" (2s)	200/sleeve, 25 sleeves/case, 5000/case
B3001-5	8-ply, 2" x 2" (5s)	200/sleeve, 25 sleeves/case, 5000/case
B3003	8-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
B3005	8-ply, 4" x 4"	200/sleeve, 20 sleeves/case, 4000/case
B3005-4	8-ply, 4" x 4" (4s)	200/sleeve, 20 sleeves/case, 4000/case
B3005-5	8-ply, 4" x 4" (5s)	200/sleeve, 20 sleeves/case, 4000/case
B3002	12-ply, 2" x 2"	200/sleeve, 40 sleeves/case, 8000/case
B3002-5	12-ply, 2" x 2" (5s)	200/sleeve, 40 sleeves/case, 8000/case
B3004	12-ply, 3" x 3"	200/sleeve, 20 sleeves/case, 4000/case
B3004-5	12-ply, 3" x 3" (5s)	200/sleeve, 20 sleeves/case, 4000/case
B3006	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case
B3006-10	12-ply, 4" x 4" (10s)	200/sleeve, 10 sleeves/case, 2000/case
B3010	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case
B3009	12-ply, 8" x 4"	200/sleeve, 10 sleeves/case, 2000/case
B3009-5	12-ply, 8" x 4" (5s)	200/sleeve, 10 sleeves/case, 2000/case
B3007-10	16-ply, 4" x 4" (10s)	200/sleeve, 10 sleeves/case, 2000/case

Sterile

B8800	8-ply, 8" x 4"	2/pouch, 25 pouches/tray, 20 trays/case, 1000/case
B2120	12-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
B3312	12-ply, 3" x 3"	1/pouch, 100 pouches/tray, 24 trays/box, 2400/case
B3120	12-ply, 3" x 3"	2/pouch, 40 pouches/tray, 30 trays/case, 2400/case
B4412	12-ply, 4" x 4"	1/pouch, 100 pouches/tray, 12 trays/box, 1200/case
B4120	12-ply, 4" x 4"	2/pouch, 25 pouches/tray, 24 trays/case, 1200/case

Continues on Page 33 >



Ritmed® Cotton Gauze Sponges Sterile

Product Code	Description	Package/Case
B4121	12-ply, 4" x 4"	10/pouch, 10 pouches/tray, 12 trays/box, 1200/case
B8120	12-ply, 8" x 4"	2/pouch, 25 pouches/tray, 20 trays/case, 1000/case
B8125	12-ply, 8" x 4"	5/pouch, 10 pouches/tray, 20 trays/case, 1000/case
B8121	12-ply, 8" x 4"	10/pouch, 10 pouches/tray, 10 trays/box, 1000/case
B2800	8-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
B4800	8-ply, 4" x 4"	2/pouch, 50 pouches/tray, 12 trays/case, 1200/case
B4181	8-ply, 4" x 4"	10/pouch, 10 pouches/tray, 12 trays/box, 1200/case



Ritmed® DisTech® Cotton Gauze Sponges Non-sterile

Product Code	Description	Package/Case
D-3001	8-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
D-3005	8-ply, 4" x 4"	200/sleeve, 20 sleeves/case, 4000/case
D-3002	12-ply, 2" x 2"	200/sleeve, 40 sleeves/case, 8000/case
D-3006	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case

Sterile

D-2800	8-ply, 2" x 2"	2/pouch, 50 pouches/tray, 30 trays/case, 3000/case
D-4800	8-ply, 4" x 4"	2/pouch, 25 pouches/tray, 24 trays/case, 1200/case
D-2212	12-ply, 2" x 2"	1/pouch, 100 pouches/tray, 24 trays/box, 2400/case
D-4412	12-ply, 4" x 4"	1/pouch, 100 pouches/tray, 12 trays/box, 1200/case
D-4120	12-ply, 4" x 4"	2/pouch, 25 pouches/tray, 24 trays/case, 1200/case



SafeGauze® Cotton Premium Non-Woven Sponges with Diamond Aperture™ Non-sterile

Product Code	Description	Package/Case
4822	4-ply, 2" x 2"	200/sleeve, 20 sleeves/case, 4000/case
4844	4-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case



SafeBasics™ All Gauze Cotton Gauze Sponges Non-sterile

Product Code	Description	Package/Case
3001	8-ply, 2" x 2"	200/sleeve, 25 sleeves/case, 5000/case
3005	8-ply, 4" x 4"	200/sleeve, 20 sleeves/case, 4000/case
3006	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case



Ritmed® Premium Cotton Gauze Sponges

X-ray detectable premium grade cotton

- Highly absorbent and low-linting.
- Latex free providing optimal patient safety and comfort.



Non-sterile, X-ray Detectable

Product Code	Description	Package/Case
6000	12-ply, 2" x 2"	200/sleeve, 40 sleeves/case, 8000/case
6003	12-ply, 3" x 3"	200/sleeve, 10 sleeves/case, 2000/case
6004	12-ply, 4" x 4"	200/sleeve, 10 sleeves/case, 2000/case
6014	12-ply, 4" x 4", Banded in 10's	200/sleeve, 10 sleeves/case, 2000/case
6008	12-ply, 8" x 4"	200/sleeve, 10 sleeves/case, 2000/case
6009	12-ply, 8" x 4", Banded in 10's	200/sleeve, 10 sleeves/case, 2000/case
6008-5	12-ply, 8" x 4", Indexed 5's	200/sleeve, 10 sleeves/case, 2000/case
6006	16-ply, 4" x 4", Banded in 10's	200/sleeve, 10 sleeves/case, 2000/case
6010	16-ply, 8" x 4"	200/sleeve, 10 sleeves/case, 2000/case
6010-10	16-ply, 8" x 4", Banded in 10's	200/sleeve, 10 sleeves/case, 2000/case

Sterile, X-ray Detectable

6128	6-ply, 8" x 4"	10/pouch, 100 pouches/case, 1000/case
6104	12-ply, 4" x 4"	10/pouch, 100 pouches/case, 1000/case
6108	12-ply, 8" x 4"	10/pouch, 100 pouches/case, 1000/case
6105	16-ply, 4" x 4"	10/pouch, 100 pouches/case, 1000/case



Ritmed® Fluff Sponges Sterile

Product Code	Description	Package/Case
A66750	4-ply, 6" x 6 ¾"	10/pouch, 48 pouches/case, 480/case
A66752	12-ply, 6" x 6 ¾"	2/pouch, 20 pouches/tray, 12 trays/case, 480/case
A66755	12-ply, 6" x 6 ¾"	5/pouch, 120 pouches/case, 600/case

 **Ritmed® Peanut Sponges in Holder**
Sterile, X-ray Detectable

Product Code	Description	Package/Case
74712	Tail, Small, $\frac{3}{8}$ "	5/pouch, 20 pouches/case, 100/case
74712X	No Tail, Small, $\frac{3}{8}$ "	5/pouch, 20 pouches/case, 100/case

 **Ritmed® Tonsil Sponges**
Sterile, X-ray Detectable, Strung on Card

Product Code	Description	Package/Case
78912	Small/Medium, $\frac{7}{8}$ "	5/pouch, 20 pouches/case, 100/case
78914	Medium, 1"	5/pouch, 20 pouches/case, 100/case

 **Ritmed® Neurological Sponges**
Non-sterile, X-ray Detectable, Single Strung

Product Code	Description	Package/Case
55992C	$\frac{1}{2}$ " x $\frac{1}{2}$ "	10/pouch, 10 pouches/case, 100/case
55993C	$\frac{1}{2}$ " x 1"	10/pouch, 10 pouches/case, 100/case
55994C	$\frac{1}{2}$ " x 3"	10/pouch, 10 pouches/case, 100/case
55995C	$\frac{3}{4}$ " x $\frac{3}{4}$ "	10/pouch, 10 pouches/case, 100/case
55996C	1" x 1"	10/pouch, 10 pouches/case, 100/case
55997C	1" x 3"	10/pouch, 10 pouches/case, 100/case

Sterile, X-ray Detectable, Single Strung

85992C	$\frac{1}{2}$ " x $\frac{1}{2}$ "	10/pouch, 10 pouches/case, 100/case
85993C	$\frac{1}{2}$ " x 1"	10/pouch, 10 pouches/case, 100/case
85994C	$\frac{1}{2}$ " x 3"	10/pouch, 10 pouches/case, 100/case
85995C	$\frac{3}{4}$ " x $\frac{3}{4}$ "	10/pouch, 10 pouches/case, 100/case
85996C	1" x 1"	10/pouch, 10 pouches/case, 100/case
85997C	1" x 3"	10/pouch, 10 pouches/case, 100/case
85998C	$\frac{1}{2}$ " x 1 $\frac{1}{2}$ "	10/pouch, 10 pouches/case, 100/case



Ritmed® Laparotomy Sponges
 Non-sterile, X-ray Dectectable, Not Washed, Tape Loop

Product Code	Description	Package/Case
1524	4-ply, 4" x 18"	Banded 5's, 200/case
1504	4-ply, 12" x 12"	Banded 5's, 200/case
1514	4-ply, 18" x 18"	Banded 5's, 200/case

Non-sterile, X-ray Dectectable, Pre-washed, Tape Loop

1501	4-ply, 12"x 12"	Banded 5's, 1400/case
1511	4-ply, 18" x 18"	Banded 5's, 800/case
1513	4-ply, 18" x 18"	Banded 5's, 200/case

Sterile, X-ray Dectectable, Not Washed, Tape Loop

1526	4-ply, 4" x 18"	5/pouch, 40 pouches/case, 200/case
1506	4-ply, 12" x 12"	5/pouch, 40 pouches/case, 200/case
1516	4-ply, 18" x 18"	5/pouch, 40 pouches/case, 200/case

Sterile, X-ray Dectectable, Pre-washed, Tape Loop

1525	4-ply, 4" x 18"	5/pouch, 40 pouches/case, 200/case
1505	4-ply, 12" x 12"	5/pouch, 40 pouches/case, 200/case
1515	4-ply, 18" x 18"	5/pouch, 40 pouches/case, 200/case

Sterile, X-ray Dectectable, Pre-washed, No Loop

1505-NL	4-ply, 12" x 12"	5/pouch, 40 pouches/case, 200/case
1515NL	4-ply, 18" x 18"	5/pouch, 40 pouches/case, 200/case



Ritmed® Conforming Bandage

Protecting wounds from further injury and contamination.

- Fits all parts of the body without restricting movement.
- Made of soft, flexible, loosely woven material.

Non-sterile

Product Code	Description	Package/Case
A205	2" x 4.1 yd	12/box, 8 boxes/case, 96/case
A305	3" x 4.1 yd	12/box, 8 boxes/case, 96/case
A405	4" x 4.1 yd	12/box, 8 boxes/case, 96/case
A605	6" x 4.1 yd	6/box, 8 boxes/case, 48/case
A405-IND	Individually Wrapped, 4" x 4.1 yd	Bulk 500/case
A605-IND	Individually Wrapped, 6" x 4.1 yd	Bulk 500/case

Sterile

A206	2" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case
A306	3" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case
A406	4" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case
A606	6" x 4.1 yd	1/pouch, 12 pouches/box, 4 boxes/case, 48/case

Ritmed® DisTech® Conforming Bandages Non-sterile

Product Code	Description	Package/Case
D205	2" x 4.1 yd	12/box, 8 boxes/case, 96/case
D305	3" x 4.1 yd	12/box, 8 boxes/case, 96/case
D405	4" x 4.1 yd	12/box, 8 boxes/case, 96/case
D605	6" x 4.1 yd	6/box, 8 boxes/case, 48/case

Sterile

D206	2" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case
D306	3" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case
D406	4" x 4.1 yd	1/pouch, 12 pouches/box, 8 boxes/case, 96/case



Product Code	Description		Package/Case
A1061-A	1" x 5 yd, Assorted		72/box
A1061-B	1" x 5 yd, Blue		72/box
A1061-T	1" x 5 yd, Tan		72/box
A2061-A	2" x 5 yd, Assorted		36/box
A2061-B	2" x 5 yd, Blue		36/box
A2061-R	2" x 5 yd, Red		36/box
A2061-T	2" x 5 yd, Tan		36/box
A2061-MIX	2" x 5 yd, Assorted		36/box
A2061-PPA	2" x 5 yd, Assorted		36/box
A3061-A	3" x 5 yd, Assorted		24/box
A3061-T	3" x 5 yd, Tan		24/box
A4061-A	4" x 5 yd, Assorted		18/box
A4061-B	4" x 5 yd, Blue		18/box
A4061-R	4" x 5 yd, Red		18/box
A4061-T	4" x 5 yd, Tan		18/box
A4061-PPA	4" x 5 yd, Assorted		18/box
A6061-A	6" x 5 yd, Assorted		12/box
A6061-T	6" x 5 yd, Tan		12/box

Not Made with Natural Rubber Latex

A1061-TLF	1" x 5 yd, Tan		72/box
A2061-ALF	2" x 5 yd, Assorted		36/box
A2061-TLF	2" x 5 yd, Tan		36/box
A3061-ALF	3" x 5 yd, Assorted		24/box
A3061-TLF	3" x 5 yd, Tan		24/box
A4061-ALF	4" x 5 yd, Assorted		18/box
A4061-TLF	4" x 5 yd, Tan		18/box
A4061-TLFB	4" x 5 yd, Tan		Bulk 100/case
A6061-TLF	6" x 5 yd, Tan		12/box
A6061-TLFB	6" x 5 yd, Tan		Bulk 72/case



Ritmed® Adhesive Bandages

Sterile

Product Code	Description	Package/Case
AF0334	Fabric, ¾" x 3", Strip	100/box
AF0310	Fabric, 1" x 3", Strip	100/box
AFF0217	Fabric, Fingertip, 1 ¾" x 2"	100/box
AFK0315	Fabric, Knuckle, 1 ½" x 3"	100/box
AF0402	Jumbo Fabric, 2" x 4", Strip	50/box
AP0402	Jumbo Plastic, 2" x 4", Strip	50/box
AP0334	Plastic, ¾" x 3", Strip	100/box
AP0310	Plastic, 1" x 3", Strip	100/box
AS0334	Sheer, ¾" x 3", Strip	100/box



Ritmed® Tubular Mesh Bandage Rolls

Not Made with Natural Rubber Latex

Product Code	Description	Package/Case
51410-LF	25 yd, Size 0, Small Finger/Toe	1 roll/box
51411-LF	25 yd, Size 1, Finger/Toe/Wrist	1 roll/box
51412-LF	25 yd, Size 2, Small Hand/Arm/Leg/Foot	1 roll/box
51413-LF	25 yd, Size 3, Medium Hand/Arm/Leg/Foot	1 roll/box
51414-LF	25 yd, Size 4, Large Hand/Arm/Leg/Foot	1 roll/box
51415-LF	25 yd, Size 5, Small Head/Shoulder/Thigh	1 roll/box
51416-LF	25 yd, Size 6, Large Head/Shoulder/Thigh	1 roll/box
51417-LF	25 yd, Size 7, Small Chest/Back/Axilla	1 roll/box
51418-LF	25 yd, Size 8, Medium Chest/Back/Axilla	1 roll/box
51419-LF	25 yd, Size 9, Large Chest/Back/Axilla	1 roll/box
51440-LF	25 yd, Size 10, X-Large Chest/Back/Axilla	1 roll/box



Ritmed® Elastic Bandage Rolls
Non-sterile, Not Made with Natural Rubber Latex

Product Code	Description	Package/Case
A620	2" x 5 yd	10/box, 5 boxes/case, 50/case
A621	3" x 5 yd	10/box, 5 boxes/case, 50/case
A622	4" x 5 yd	10/box, 5 boxes/case, 50/case
A623	6" x 5 yd	10/box, 5 boxes/case, 50/case

Non-sterile, Not Made with Natural Rubber Latex, Velcro

A632	2" x 5 yd	10/box, 5 boxes/case, 50/case
A633	3" x 5 yd	10/box, 5 boxes/case, 50/case
A634	4" x 5 yd	10/box, 5 boxes/case, 50/case
A636	6" x 5 yd	10/box, 5 boxes/case, 50/case



SafeBasics™ Triangular Bandages
Non-sterile

Product Code	Description	Package/Case
5460	40" x 40" x 56"	1/pouch, 12 pouches/bag, 20bags/case, 240/case



Ritmed® Abdominal Pads

*Provides cushioning and protection without adhering to wound.
[Sterile and non-sterile]*

- Highly absorbent, soft material construction for patient comfort.
- Disperse fluids quickly and evenly and sealed on all 4 sides.



Non-sterile

Product Code	Description	Package/Case
A7159	5" x 9"	576/bag, 1 bag/case, 576/case
A7181	8" x 10"	Bulk, 432/bag, 1 bag/case, 432/case
A7184	8" x 24"	144/bag, 1 bag/case, 144/case
A7182	8" x 20 yd	1 roll/bag, 6 bags/case, 6/case
Sterile		
A7059	5" x 9"	1/pouch, 20 pouches/tray, 20 trays/case, 400/case
A7087	8" x 7 1/2"	1/pouch, 12 pouches/tray, 20 trays/case, 240/case
A7081	8" x 10"	1/pouch, 24 pouches/tray, 15 trays/case, 360/case
A7084	8" x 24"	1/pouch, 10 pouches/tray, 6 trays/case, 60/case
A7011	10" x 12"	1/pouch, 10 pouches/tray, 20 trays/case, 200/case



Ritmed® Non-Adherent Pads
Sterile

Product Code	Description	Package/Case
A1923	2" x 3"	1/pouch, 100 pouches/box, 12 boxes/case, 1200/case
A1934	3" x 4"	1/pouch, 100 pouches/box, 12 boxes/case, 1200/case
A1938	3 1/2" x 8"	1/pouch, 50 pouches/box, 12 boxes/case, 600/case



Ritmed® Disposable Underpad

Product Code	Description	Package/Case
A7200	3-ply, 12" x 17"	Bulk 500/case
A7201	3-ply, 12" x 17"	50/box, 20 boxes/case, 1000/case
A7202	3-ply, 17" x 24"	Bulk 300/case
A7203	3-ply, 17" x 24"	30/box, 10 boxes/case, 300/case



Ritmed® Eye Pad Dressing Non-Woven Sponges Sterile

Product Code	Description	Package/Case
A1110	1 5/8" x 2 5/8"	1/pouch, 50 pouches/tray, 12 trays/case, 600/case
A1111	2 1/8" x 2 5/8"	1/pouch, 50 pouches/tray, 12 trays/case, 600/case



Ritmed® Surgical Tape Non-sterile

Product Code	Description	Package/Case
A5010	Clear, 1" x 10 yd	12/box, 12 boxes/case, 144/case
A5010-15FT	Clear, 1" x 15'	Bulk 144/case
A5020	Clear, 2" x 10 yd	6/box, 12 boxes/case, 72/case
A5105	Paper, 1/2" x 10 yd	24/box, 12 boxes/case, 288/case
A5110	Paper, 1" x 10 yd	12/box, 12 boxes/case, 144/case
A5110-15FT	Paper, 1" x 15'	144/case
A5120	Paper, 2" x 10 yd	6/box, 12 boxes/case, 72/case
A5205	Silk, 1/2" x 10 yd	24/box, 12 boxes/case, 288/case
A5210	Silk, 1" x 10 yd	12/box, 12 boxes/case, 144/case
A5220	Silk, 2" x 10 yd	6/box, 12 boxes/case, 72/case
A5230	Silk, 3" x 10 yd	4/box, 12 boxes/case, 48/case



Ritmed® Vital Gauze Bandage Rolls Non-sterile

Product Code	Description	Package/Case
A6903	2-ply, 3" x 131"	12/bag, 8 bags/case, 96/case
A6904	2-ply, 4" x 131"	12/bag, 8 bags/case, 96/case

Sterile

A6924	2-ply, 4" x 131"	1/pouch, 12 pouches/bag, 8 bags/case, 96/case
-------	------------------	---



Ritmed® Krimped Gauze Rolls Non-sterile

Product Code	Description	Package/Case
A454	6-ply, 4 1/2" x 4.1 yd	100/case

Sterile

A455	Premium, 6-ply, 4 1/2" x 4.1 yd	1/pouch, 100 pouches/case, 100/case
D455	Distech, 6-ply, 4 1/2" x 4.1 yd	1/pouch, 100 pouches/case, 100/case



Ritmed® Transparent Dressing Sterile

Product Code	Description	Package/Case
A60044	With window and label to write on, 4" x 4 ¾"	1/pouch, 50 pouches/box, 4 boxes/case, 200/case



Ritmed® AssureGauze Packing Gauze Sterile, X-ray Detectable, Non-strung

Product Code	Description	Package/Case
84442	2" x 36"	1/pouch, 100 pouches/case, 100/case



Ritmed® Kittner Dissector Sponges Sterile, X-ray Detectable

Product Code	Description	Package/Case
73512	9/16" x 1/4"	5/pouch, 100/case



SafeGauze® HemoStat™ Topical Hemostatic Dressing Sterile

Product Code	Description	Package/Case
4900-STATC	4-ply, ¾" x ¾"	20/box



Countway Systems

Product Code	Description	Package/Case
2212A	Round Countway Stand	Wire basket
2212B	Square Countway Stand	Wire basket
2213	5 Pocket Sponge Countway Panel, 29" x 14 ½" Green	50/box, 8 boxes/case, 400/case
2213-BL	5 Pocket Sponge Countway Panel, 29" x 14 ½" Blue	50/box, 8 boxes/case, 400/case
2214	Inner Center Sponge Collecting Bag	50/box, 8 boxes/case, 400/case



Sterilization

Peace of mind.

Our reliable self-sealing sterilization pouches, CSR sterilization wraps, and sterilization reels provide the confidence that comes from knowing that your instruments have been thoroughly sterilized and are safe to reuse.



SafeBasics™ Self-Sealing Sterilization Pouches

Product Code	Description	Package/Case
99000	2 ¾" x 9"	200/box, 6 boxes/case, 1200/case
99005	3 ½" x 5¼"	200/box, 5 boxes/case, 1000/case
99010	3 ½" x 9"	200/box, 6 boxes/case, 1200/case
99015	2 ¼" x 4"	200/box, 10/boxes/case, 2000/case
99025	5 ¼" x 10"	200/box, 10 boxes/case, 2000/case
99030	7 ½" x 13"	200/box, 5 boxes/case, 1000/case
99010-B	3 ½" x 9"	Bulk 500/box, 6 boxes/case, 3000/case



SafeSeal® Quattro Self-Sealing Sterilization Pouches with TruePress™ Technology

Product Code	Description	Package/Case
88000-4	2 ¾" x 9"	200/box, 6 boxes/case, 1200/case
88005-4	3 ½" x 5 ¼"	200/box, 5 boxes/case, 1000/case
88010-4	3 ½" x 9"	200/boxes, 6 boxes/case, 1200 case
88010-4B	3 ½" x 9"	Bulk 500/box, 6 boxes/case, 3000 case
88015-4	2 ¼" x 4"	200/box, 10 boxes/case, 2000 case
88020-4	5 ¼" x 6 ½"	200/box, 10 boxes/case, 2000 case
88025-4	5 ¼" x 10"	200/box, 10boxes/case, 2000 case
88030-4	7 ½" x 13"	200/box, 5 boxes/case, 1000 case
88035-4	10" x 14"	200/box, 5 boxes/case, 1000 case
88040-4	12" x 17"	200/box, 5 boxes/case, 1000 case
88045-4	4 ¼" x 11"	200/box, 6 boxes/case, 1200 case



Peelable Sterilization Reels

Product Code	Description	Package/Case
9510	2" x 100'	10 rolls/case
9520	3" x 100'	10 rolls/case
9530	4" x 100'	10 rolls/case
9540	6" x 100'	10 rolls/case




Hedy® Quadra Fold SMMS Wraps

Product Code	Description	Package/Case
620WR-15MW	15" x 15"	250/pack, 1 pack/case, 250/case
620WR-20MW	20" x 20"	125/pack, 1 pack/case, 125/case
620WR-24MW	24" x 24"	125/pack, 1 pack/case, 125/case



Hedy® Steam Autoclave Tapes

Product Code	Description	Package/Case
630TA-1	1" x 55 yd	1/case



Essential Supplies

Little things that go a long way.

Whether swabs, applicators, tongue depressors, or dry or wet wipes, these essentials are a must for your daily needs.



SafeBasics™ Cotton Tipped Applicators

Non-sterile

Product Code	Description	Package/Case
803	3"	100/bag, 10 bags/box, 10 boxes/case, 10000/case
806	6"	100/bag, 10 bags/box, 10 boxes/case, 10000/case

Sterile

804	3"	2/pouch, 100 pouches/box, 10 boxes/case, 2000/case
807	6"	2/pouch, 100 pouches/box, 10 boxes/case, 2000/case



Solon™ Swabs & Applicators

Non-sterile

Product Code	Description	Package/Case
56101	Cotton Tip Applicator Plastic Stem, 3"	100/bag, 10 bags/box, 10 boxes/case, 10000/case
56122	Cotton Tip Applicator Tapered Wood Stem, 3"	100/bag, 10 bags/box, 10 boxes/case, 10000/case
56201	Cotton Tip Applicator Plastic Stem, 6"	100/bag, 10 bags/box, 10 boxes/case, 10000/case
56225	Cotton Tip Applicator Tapered Wood Stem, 6"	100/bag, 10 bags/box, 10 boxes/case, 10000/case
56701	Cotton Tip Applicator Plastic Stem, 6"	1/pouch, 100 pouches/box, 10 boxes/case, 1000/case
57500	Proctoscopic Applicator, Rayon Tip, Plastic Stem, 16"	50/box, 10 boxes/case, 500/case
57600	OB/GYN Applicator Rayon Tip Plastic Stem, 8"	50/box, 10 boxes/case, 500/case

Sterile

4001D	Polyester Tip, Plastic Stem, 6"	1/pouch, 100 pouches/box, 10 boxes/case, 1000/case
4001R	Rayon Tip, Plastic Stem, 6"	1/pouch, 100 pouches/box, 10 boxes/case, 1000/case
4002D	Polyester Tip, Plastic Stem, 6"	2/pouch, 100 pouches/box, 10 boxes/case, 2000/case
56801	Cotton Tip Applicator, Plastic Stem, 6"	2/pouch, 100 pouches/box, 10 boxes/case, 2000/case
57603	OB/GYN Applicator Rayon Tip, Paper Stem, 8"	2/pouch, 50 pouches/box, 10 boxes/case, 1000/case
57604	Obgyn Swab Rayon Tip, Plastic Stem, 8"	Bulk 500/case



Paragon® Disposable Scalpel, Sterile

Product Code	Description	Package/Case
90000-10	#10	1/pouch, 10 pouches/case, 10/case
90000-11	#11	1/pouch, 10 pouches/case, 10/case
90000-12	#12	1/pouch, 10 pouches/case, 10/case
90000-15	#15	1/pouch, 10 pouches/case, 10/case
90000-20	#20	1/pouch, 10 pouches/case, 10/case
90000-21	#21	1/pouch, 10 pouches/case, 10/case
90000-22	#22	1/pouch, 10 pouches/case, 10/case
90000-23	#23	1/pouch, 10 pouches/case, 10/case
90000-24	#24	1/pouch, 10 pouches/case, 10/case



Paragon® Disposable Blades, Sterile Stainless Steel

Product Code	Description	Package/Case
90010-10	#10	1/pouch, 100 pouches/case, 100/case
90010-11	#11	1/pouch, 100 pouches/case, 100/case
90010-12	#12	1/pouch, 100 pouches/case, 100/case
90010-15	#15	1/pouch, 100 pouches/case, 100/case
90010-20	#20	1/pouch, 100 pouches/case, 100/case
90010-21	#21	1/pouch, 100 pouches/case, 100/case
90010-22	#22	1/pouch, 100 pouches/case, 100/case
90010-23	#23	1/pouch, 100 pouches/case, 100/case
90010-24	#24	1/pouch, 100 pouches/case, 100/case

Carbon Steel

90011-10	#10	1/pouch, 100 pouches/case, 100/case
90011-11	#11	1/pouch, 100 pouches/case, 100/case
90011-12	#12	1/pouch, 100 pouches/case, 100/case
90011-15	#15	1/pouch, 100 pouches/case, 100/case
90011-20	#20	1/pouch, 100 pouches/case, 100/case
90011-21	#21	1/pouch, 100 pouches/case, 100/case
90011-22	#22	1/pouch, 100 pouches/case, 100/case
90011-23	#23	1/pouch, 100 pouches/case, 100/case
90011-24	#24	1/pouch, 100 pouches/case, 100/case



O.R. Towels, Pre-washed Non-sterile

Product Code	Description	Package/Case
5100-B	17" x 25", Blue	<input checked="" type="radio"/> Bulk 100/case
5100-G	17" x 25", Green	<input checked="" type="radio"/> Bulk 100/case
5400-B	17" x 25", Blue	<input checked="" type="radio"/> Bulk 400/case
5400-W	17" x 25", White	<input type="radio"/> Bulk 400/case

Sterile

5001-B	17" x 25", Blue	<input checked="" type="radio"/> 1/pouch, 80 pouches/case, 80/case
5002-B	17" x 25", Blue	<input checked="" type="radio"/> 2/pouch, 40 pouches/case, 80/case
5002-BX	X-ray Detectable, 17" x 25", Blue	<input checked="" type="radio"/> 2/pouch, 40 pouches/case, 80/case
5002-GX	X-ray Detectable, 17" x 25", Green	<input checked="" type="radio"/> 2/pouch, 40 pouches/case, 80/case
5002-W	17" x 25", White	<input type="radio"/> 2/pouch, 40 pouches/case, 80/case
5002-WX	X-ray Detectable, 17" x 25", White	<input type="radio"/> 2/pouch, 40 pouches/case, 80/case
5004-B	17" x 25", Blue	<input checked="" type="radio"/> 4/pouch, 20 pouches/case, 80/case
5004-BX	X-ray Detectable, 17" x 25", Blue	<input checked="" type="radio"/> 4/pouch, 20 pouches/case, 80/case
5004-G	17" x 25", Green	<input checked="" type="radio"/> 4/pouch, 20 pouches/case, 80/case
5004-W	17" x 25", White	<input type="radio"/> 4/pouch, 20 pouches/case, 80/case
5004-WX	X-ray Detectable, 17" x 25", White	<input type="radio"/> 4/pouch, 20 pouches/case, 80/case
5006-B	17" x 25", Blue	<input checked="" type="radio"/> 6/pouch, 12 pouches/case, 72/case
5006-G	17" x 25", Green	<input checked="" type="radio"/> 6/pouch, 12 pouches/case, 72/case
5006-GX	X-ray Detectable, 17" x 25", Green	<input checked="" type="radio"/> 6/pouch, 12 pouches/case, 72/case
5008-B	17" x 25", Blue	<input checked="" type="radio"/> 8/pouch, 10 pouches/case, 80/case
5010-B	17" x 25", Blue	<input checked="" type="radio"/> 10 pouch, 8 pouches/case, 80/case



Professional Towel [Made in USA]

Product Code	Description	Package/Case
8190	3-ply, 13" x 18", White	<input type="radio"/> 500/case



Pre-Moistened Washcloths, Spunlace

Product Code	Description	Package/Case
A40001	Hard Tub, 9" x 13"	50/tub, 12 tubs/case, 600/case
A40002	9" x 13"	64/pack, 8 packs/case, 512/case
A40002-SF	Unscented, 9" x 13"	64/pack, 8 packs/case, 512/case
A40003	9" x 13"	50/pack, 12 packs/case, 600/case



General Purpose Drape Sheets, Non-sterile

Product Code	Description	Package/Case
A8400	36" x 72", White	○ 25/bag, 10 bag/case, 250/case



SafeBasics Dry-Back Bibs

Product Code	Description	Package/Case
8280	2-ply, Lavender	● 500/case
8281	2-ply, Silver Grey	● 500/case
8282	2-ply, Blue	● 500/case
8283	2-ply, White	○ 500/case
8284	2-ply, Yellow	● 500/case
8285	2-ply, Green	● 500/case
8286	2-ply, Dusty Rose	● 500/case
8287	2-ply, Beige	● 500/case
8288	2-ply, Peach	● 500/case
8289	2-ply, Aqua	● 500/case
8292	2-ply, Black	● 500/case
8380	3-ply, Lavender	● 500/case
8381	3-ply, Silver Grey	● 500/case
8382	3-ply, Blue	● 500/case
8383	3-ply, White	○ 500/case
8384	3-ply, Yellow	● 500/case
8385	3-ply, Green	● 500/case
8386	3-ply, Dusty Rose	● 500/case
8387	3-ply, Beige	● 500/case
8388	3-ply, Peach	● 500/case
8389	3-ply, Aqua	● 500/case



Pillow Cases

Product Code	Description	Package/Case
A1090	Disposable, 21" x 30", White	<input type="radio"/> 100/case



SafeBasics™ Examination Table Paper

Product Code	Description	Package/Case
80201	21" x 225', Smooth White	<input type="radio"/> 12 rolls/case
80202	21" x 125', Crepe White	<input type="radio"/> 12 rolls/case
80203	18" x 225', Smooth White	<input type="radio"/> 12 rolls/case
80204	18" x 125', Crepe White	<input type="radio"/> 12 rolls/case
80205	Individually Wrapped, 24" x 225', Smooth White	<input type="radio"/> 12 rolls/case



SafeBasics™ Tongue Depressors Non-sterile

Product Code	Description	Package/Case
40006	6"	500/box, 10 boxes/case, 5000/case

Sterile

40004	5 ½"	1/pouch, 100 pouches/box, 10 boxes/case, 1000/case
40007	6"	1/pouch, 100 pouches/box, 10 boxes/case, 1000/case



Tongue Depressors Non-sterile

Product Code	Description	Package/Case
50000	SolonCare, Junior, 5 ½"	500/box, 10 boxes/case, 5000/case
55900	Solon Wood Applicator, Standard, 6"	1000/box, 20 boxes/case, 20000/case



SafeBasics™ Plastic Cups

Product Code	Description	Package/Case
105	5 oz, Green	 100/bag, 10 bags/case, 1000/case
106	5 oz, Lavender	 100/bag, 10 bags/case, 1000/case
107	5 oz, Dusty Rose	 100/bag, 10 bags/case, 1000/case
108	5 oz, Yellow	 100/bag, 10 bags/case, 1000/case
110	5 oz, White	 100/bag, 10 bags/case, 1000/case
112	5 oz, Blue	 100/bag, 10 bags/case, 1000/case
121	5 oz, Aqua	 100/bag, 10 bags/case, 1000/case



SafeBasics™ Poly-Coated Paper Cups

Product Code	Description	Package/Case
114-CH	4 oz, Healthy Teeth	 100/bag, 10 bags/case, 1000/case
113-CH	4 oz, Bubbles Design	 100/bag, 10 bags/case, 1000/case
115-CH	5 oz, Bubbles Design	 100/bag, 10 bags/case, 1000/case
116-CH	5 oz, Healthy Teeth	 100/bag, 10 bags/case, 1000/case



SafePure® 80% Liquid Hand Sanitizer

Product Code	Description	Package/Case
7000	Cartridge, 1L	12 units/case
7100	Bottle, 1L	6 units/case
7500	Bottle, 500ml	12 units/case



Isopropyl Alcohol

Product Code	Description	Package/Case
30070	4L, 70%	4 bottles/case
30099	4L, 99%	4 bottles/case



Dry Wipes

Product Code	Description	Package/Case
A40010	Airlaid Premium, ¼ Fold, 10" x 13"	50/pack, 10 packs/case, 500/case
A40013	Airlaid Premium, 12" x 13"	50/pack, 10 packs/case, 500/case
A41013-1	Airlaid Premium, Soft, 10" x 13"	50/pack, 20 packs/case, 1000/case
BP-40010	Airlaid, Soft, 9" x 12"	50/pack, 10 packs/case, 500/case
BP40013	Airlaid, Soft, 12" x 13"	50/pack, 10 packs/case, 500/case
CM-40010	DRC "Cellumed" High Absorbency, 10" x 13"	50/pack, 10 packs/case, 500/case
CM-40013	DRC "Cellumed" High Absorbency, 12" x 13"	50/pack, 10 packs/case, 500/case
CM-41213-1	DRC "Cellumed" High Absorbency, ¼ Fold, 12" x 13"	50/pack, 20 packs/case, 1000/case
HC40010	Hydro-Cell Strong, Soft, 10" x 13"	50/pack, 10 packs/case, 500/case
HC41213-1	Hydro-Cell Strong, Soft, 12" x 13"	50/pack, 20 packs/case, 1000/case
FS-41013	Food Service, Rayon, 10" x 13"	50/pack, 10 packs/case, 500/case
904014	Dispersible Fluff Spun, 9" x 13"	Bulk 1000/case
SP-40813	Spunlace, 8" x 13"	50/pack, 10 packs/case, 500/case
SP-41013	Spunlace, 10" x 13"	50/pack, 10 packs/case, 500/case
SP-41213-1	Spunlace, 12" x 13"	50/pack, 20 packs/case, 1000/case
QS-41013	Spunlace Premium, 10" x 13"	50/pack, 10 packs/case, 500/case
QS-41213-1	Spunlace Premium, 12" x 13"	50/pack, 20 packs/case, 1000/case
RS-41213	Scrim Reinforced White Towel, 4-ply, 12" x 13"	50/pack, 18 packs/case, 900/case



Flushable Wet Wipes

Product Code	Description	Package/Case
A40008	Hard Tub, 9" x 13"	60/tub, 8 tubs/case, 480/case
A40028	9" x 13"	24/pack, 24 packs/case, 576/case



Ensuring healthier lives and a safer world for all.

Our responsiveness to the everchanging infection prevention needs and our commitment to reliable, high-quality products helps ensure you have the protection you need, when you need it most.



For more information, please visit medicom.com or contact us at: cservice@medicom.ca / 1-800-361-2862

AMD Medicom Inc. 2555 Chemin de l'Aviation, Montreal, Quebec H9P 2Z2 Canada

©2023 AMD Medicom Inc. All rights reserved.

AUTHORIZED DISTRIBUTOR



Rev. 00