

**WE HAVE YOU
COVERED**



Reliable protection without compromise to comfort and style

Reduce the risk of cross-contamination and eliminate the need for laundering with our portfolio of high-quality, single-use PPE options.

SafeBasics

SafeBasics protective apparel provides a basic level of protection at an attractive price.

- Combines protection and comfort
- Versatile design and material options
- Cost effective



Isolation Gowns

Single-use gown for procedures requiring protection from spray and splatter. Fluid-repellent design provides added protection against cross-contamination.



Lab Gowns

Front snap closure, knitted cuffs and collar for comfortable, easy on-and-off fit.



Fluid Impervious Gowns

Polyethylene coating provides fluid resistance for procedures involving the potential for cross-contamination.



Item	Description	Unit/Pkg	Pkg/Case
D-8009	Isolation Gown M/L - Yellow	10	5
D-8010	Isolation Gown XL - Yellow	10	5
8023	Isolation Gown Reg - Blue	10	5

Item	Description	Unit/Pkg	Pkg/Case
8016	Lab Gown S - Sky Blue	10	5
8017	Lab Gown M - Sky Blue	10	5
8018	Lab Gown L - Sky Blue	10	5
8019	Lab Gown XL - Sky Blue	10	5

Item	Description	Unit/Pkg	Pkg/Case
8022	Fluid Impervious Gown Reg - White	10	5
8024	Fluid Impervious Gown XL - White	10	5

Bouffant Caps

Hair and head fluid protection with breathable, low-linting design.

Beard Covers

Face protection with an elastic edge and headband loop for optimal fit and comfort.



Item	Description	Unit/Pkg	Pkg/Case
8000	Bouffant Caps 21" - White	100	5
8001	Bouffant Caps 21" - Blue	100	5
8003	Bouffant Caps 24" - White	100	5
8004	Bouffant Caps 24" - Blue	100	5
8035	Beard Covers - White	100	10

Shoe Covers

Fluid repellent footwear protection.



Item	Description	Unit/Pkg	Pkg/Case
8006	Shoe Cover Non-Conductive Non-Skid - XL - Blue	100	3
8007	Shoe Cover Non-Conductive Non-Skid - Reg - Blue	100	3
8008	Shoe Cover Non-Conductive Reg - Blue	100	3
8036	Shoe Cover Non-Conductive Non-Skid - XL - White	100	3

SafeWear™

Premium Protective Apparel

- Lightweight, soft, supple fabrics keep wearers cool and comfortable.
- Secure front snaps.
- Latex free.
- Made from AAMI PB70 Level 2 tested material.
- Complies with recommended OSHA best practices.
- Available in a wide range of sizes, from XS - 2XL.



Form-Fit Isolation Gown™
Premium isolation gown

Item	Description	Unit/ Pkg	Pkg/ Case
8113 ¹	Sunny Yellow	12	5
8114 ¹	Bright Blue	12	5



Hipster Jacket™
Premium hip length jacket

Item	Description	Unit/ Pkg	Pkg/ Case
8101*	Deep Blue	12	5
8103*	Plum Purple	12	5
8104	Soft Blue	12	5
8105	White Frost	12	5



High Performance Lab Coat™
Premium lab coat

Item	Description	Unit/ Pkg	Pkg/ Case
8107	Deep Blue	12	5
8109	Plum Purple	12	5
8110	White Frost	12	5
8112	Soft Blue	12	5

¹Available only in regular size.

Size:

XS: Extra Small

A: Small

B: Medium

C: Large

D: Extra Large

E: 2X Large

*Available only in Hipster Jacket, Deep blue & Plum Purple.

Hipster Jackets and High Performance Lab Coats: Small, Medium, Large, and Extra Large available in all colours. 2X Large available in Deep Blue (only).

Protective Apparel Best Practices

DONNING

- Perform hand hygiene before removing item from package.
- Select appropriate garment for intended use (length, coverage, fluid resistance, etc.)
- Fasten correctly.

DOFFING

- Remove apparel immediately after use and dispose in a designated receptacle. Do not re-use.
- Perform proper hand hygiene after removing apparel to prevent cross-contamination.

KEEP IN MIND:

- Protective apparel should not be worn in restrooms, personal offices or break rooms.
- Change protective apparel between patients and as soon as possible when penetrated by blood or other potentially infectious fluids.
- Turn garment inside out when removing to keep infectious material inside.



[medicom.com](https://www.medicom.com) | 1-800-361-2862

Rev.02
©2024 AMD Medicom Inc. All rights reserved.



AUTHORIZED DISTRIBUTOR

